




PTPN22 Antibody

Product Code	CSB-PA780242
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2R2
Immunogen	Synthesized peptide derived from Internal of Human PTPN22.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PTPN22 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Hematopoietic cell protein-tyrosine phosphatase 70Z-PEP; Lymphoid phosphatase; lymphoid-specific protein tyrosine phosphatase; LyP; LYP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTPN22

Image



Western blot analysis of extracts from HeLa/Jurkat cells, using PTPN22 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.