




USP11 Antibody

Product Code	CSB-PA986850
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51784
Immunogen	Synthesized peptide derived from Internal of Human USP11.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total USP11 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Deubiquitinating enzyme 11; Ubiquitin carboxyl-terminal hydrolase 11; ubiquitin carboxyl-terminal hydrolase; X-linked; ubiquitin specific peptidase 11
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP11

Image



Western blot analysis of extracts from 293 cells, using USP11 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.