




USP15 Antibody

Product Code	CSB-PA040632
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4E8
Immunogen	Synthesized peptide derived from N-terminal of Human USP15.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total USP15 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 15; Ubiquitin thioesterase 15; Ubiquitin-specific-processing protease 15; Deubiquitinating enzyme 15; Unph-2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP15

Image



Western blot analysis of extracts from HeLa cells, using USP15 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.