





USP19 Antibody

Product Code	CSB-PA173374
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94966
Immunogen	Synthesized peptide derived from internal of Human USP19.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total USP19 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	deubiquitinating enzyme 19; EC 3.1.2.15; KIAA0891; UB19; ubiquitin carboxyl-terminal hydrolase 19
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP19

Image



Western blot analysis of extracts from 3T3 cells, using USP19 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using USP19 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.