




USP30 Antibody

Product Code	CSB-PA067099
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70CQ3
Immunogen	Synthesized peptide derived from internal of Human USP30.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total USP30 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	flj40511; mgc10702; ubiquitin specific peptidase 30; ubiquitin-specific proteinase 30; ubp30
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP30

Image



Western blot analysis of extracts from HeLa cells, using USP30 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.