




USP38 Antibody

Product Code	CSB-PA578250
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NB14
Immunogen	Synthesized peptide derived from internal of Human USP38.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total USP38 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 38; Ubiquitin thioesterase 38; Ubiquitin-specific-processing protease 38; Deubiquitinating enzyme 38; HP43.8KD
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP38

Image



Western blot analysis of extracts from 293 cells, using USP38 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.