




USP42 Antibody

Product Code	CSB-PA151486
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9J4
Immunogen	Synthesized peptide derived from internal of Human USP42.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total USP42 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	deubiquitinating enzyme 42; EC 3.1.2.15; UB42; ubiquitin carboxyl-terminal hydrolase 42; ubiquitin thiolesterase 42
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP42

Image



Western blot analysis of extracts from HT-29 cells, using USP42 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.