




USP45 Antibody

Product Code	CSB-PA226027
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70EL2
Immunogen	Synthesized peptide derived from internal of Human USP45.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total USP45 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 45; Ubiquitin thioesterase 45; Ubiquitin-specific-processing protease 45; Deubiquitinating enzyme 45; USP45
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP45

Image



Western blot analysis of extracts from A549 cells, using USP45 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.