




UBAC1 Antibody

Product Code	CSB-PA958011
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSL1
Immunogen	Synthesized peptide derived from internal of Human UBAC1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total UBAC1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase KPC2; glialblastoma cell differentiation-related protein 1; Kip1 ubiquitination-promoting complex 2; UBAC1; UBADC1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBAC1

Image



Western blot analysis of extracts from 293 cells and HUVEC cells, using UBAC1 antibody.