




UBAP2L Antibody

Product Code	CSB-PA261751
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14157
Immunogen	Synthesized peptide derived from internal of Human UBAP2L.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total UBAP2L protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	UBAP2L; ubiquitin-associated protein 2-like; UBP2L;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBAP2L

Image



Western blot analysis of extracts from HT-29 cells and 293 cells, using UBAP2L antibody.