




UBE2D2 Antibody

Product Code	CSB-PA826958
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62837
Immunogen	Synthesized peptide derived from C-terminal of Human UBE2D2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total UBE2D2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 D2; Ubiquitin carrier protein D2; Ubiquitin-conjugating enzyme E2(17)KB 2; Ubiquitin-conjugating enzyme E2-17 kDa 2; Ubiquitin-protein ligase D2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBE2D2

Image



Western blot analysis of extracts from HeLa cells, A549 cells and K562 cells, using UBE2D2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.