




UBE2T Antibody

Product Code	CSB-PA212556
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPD8
Immunogen	Synthesized peptide derived from internal of Human UBE2T.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total UBE2T protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 T; Ubiquitin-protein ligase T; Ubiquitin carrier protein T; UBE2T;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBE2T

Image



Western blot analysis of extracts from HeLa cells, Jurkat cells and K562 cells, using UBE2T antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.