





THOC5 Antibody

Product Code	CSB-PA890051
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13769
Immunogen	Synthesized peptide derived from internal of Human THOC5.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total THOC5 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	THO complex subunit 5 homolog; NF2/meningioma region protein pK1.3; Placental protein 39.2; THOC5; C22orf19
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	THOC5

Image



Western blot analysis of extracts from JK cells and HepG2 cells, using THOC5 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using THOC5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.