




USP6NL Antibody

Product Code	CSB-PA910583
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92738
Immunogen	Synthesized peptide derived from internal of Human USP6NL.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total USP6NL protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	USP6 N-terminal-like protein; Related to the N-terminus of tre; RN-tre; USP6NL; KIAA0019
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP6NL

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using USP6NL antibody.