




ATP6V1C2 Antibody

Product Code	CSB-PA985340
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEY4
Immunogen	Synthesized peptide derived from internal of Human ATP6V1C2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ATP6V1C2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Vacuolar proton pump subunit C 2; V-ATPase subunit C 2; ATP6V1C2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP6V1C2

Image



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using ATP6V1C2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.