




VTI1B Antibody

Product Code	CSB-PA153198
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UEU0
Immunogen	Synthesized peptide derived from N-terminal of Human VTI1B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total VTI1B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	v-SNARE; Vesicle transport through interaction with t-SNAREs homolog 1B; vesicle transport through interaction with t-SNAREs homolog 1B (yeast); Vesicle transport v-SNARE protein Vti1-like 1; vesicle-associated soluble NSF attachment protein receptor
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VTI1B

Image



Western blot analysis of extracts from HuvEc cells, using VTI1B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.