





ZBTB40 Antibody

Product Code	CSB-PA030777
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUA8
Immunogen	Synthesized peptide derived from C-terminal of Human ZBTB40.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZBTB40 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger and BTB domain-containing protein 40; ZBTB40; KIAA0478;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZBTB40

Image



Western blot analysis of extracts from 293 cells, using ZBTB40 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using ZBTB40 antibody.