




ZC3H11A Antibody

Product Code	CSB-PA782142
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75152
Immunogen	Synthesized peptide derived from internal of Human ZC11A.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ZC11A protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ZC3H11A; ZC11A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H11A

Image



Western blot analysis of extracts from COS cells, using ZC11A antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.