




ZC3H13 Antibody

Product Code	CSB-PA921508
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T200
Immunogen	Synthesized peptide derived from internal of Human ZC3H13.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ZC3H13 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ZC313;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H13

Image



Western blot analysis of extracts from HepG2 cells, A549 cells and Jurkat cells, using ZC3H13 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.