





ZC3H8 Antibody

Product Code	CSB-PA832766
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5P1
Immunogen	Synthesized peptide derived from internal of Human ZC3H8.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZC3H8 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger CCCH domain-containing protein 8; ZC3H8; ZC3HDC8;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H8

Image



Western blot analysis of extracts from HepG2 cells, using ZC3H8 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using ZC3H8 antibody.