




CLNS1A Antibody

Product Code	CSB-PA159505
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54105
Immunogen	Synthesized peptide derived from C-terminal of Human CLNS1A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CLNS1A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Methylosome subunit pICln; Chloride conductance regulatory protein ICln; I(Cln); Chloride channel; nucleotide sensitive 1A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLNS1A

Image



Western blot analysis of extracts from HuvEc cells, K562 cells, HeLa cells, HepG2 cells, A549 cells and Jurkat cells, using CLNS1A antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.