




PTGR2 Antibody

Product Code	CSB-PA234635
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8N7
Immunogen	Synthesized peptide derived from internal of Human ZADH1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ZADH1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Prostaglandin reductase 2; PRG-2; EC 1.3.1.48; 15-oxoprostaglandin 13-reductase; Zinc-binding alcohol dehydrogenase domain-containing protein 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTGR2

Image



Western blot analysis of extracts from HUVEC cells, K562 cells and HT-29 cells, using ZADH1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.