




ZNHIT6 Antibody

Product Code	CSB-PA259510
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWK9
Immunogen	Synthesized peptide derived from internal of Human CA181.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CA181 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HIT-type zinc finger-containing protein C1orf181; HIT-type zinc finger-containing protein C1orf181; C1orf181;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNHIT6

Image



Western blot analysis of extracts from Jurkat cells and HeLa cells, using CA181 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.