




TBX3 Antibody

Product Code	CSB-PA052231
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15119
Immunogen	Synthesized peptide derived from internal of Human TBX3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TBX3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	T-box protein 3; T-box transcription factor TBX3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TBX3

Image



Western blot analysis of extracts from A549 cells, using TBX3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.