





ADCY9 Antibody

Product Code	CSB-PA178331
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60503
Immunogen	Synthesized peptide derived from N-terminal of Human ADCY9.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ADCY9 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP pyrophosphate-lyase 9; Adenylate cyclase type IX; Adenylyl cyclase 9; mKIAA0520;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADCY9

Image



Western blot analysis of extracts from HuvEc cells, using ADCY9 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using ADCY9 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.