





TBPL2 Antibody

Product Code	CSB-PA264692
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6SJ96
Immunogen	Synthesized peptide derived from internal of Human TBPL2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TBPL2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TATA box-binding protein-like protein 2; TBP-like protein 2; TATA box-binding protein-related factor 3; TBP-related factor 3; TBPL2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TBPL2

Image



Western blot analysis of extracts from HeLa cells, using TBPL2 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TBPL2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.