




LIPB1 Antibody

Product Code	CSB-PA222467
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C20615
Immunogen	Synthesized peptide derived from internal of Human LIPB1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total LIPB1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Liprin-beta-1; Protein tyrosine phosphatase receptor type f polypeptide-interacting protein-binding protein 1; PTPRF-interacting protein-binding protein 1; hSGT2; PPFIBP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LIPB1

Image



Western blot analysis of extracts from A549 cells, Jurkat cells, HeLa cells, HepG2 cells, COLO205 cells and HuvEc cells, using LIPB1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.