




RNF113B Antibody

Product Code	CSB-PA576681
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZP6
Immunogen	Synthesized peptide derived from Internal of Human RNF113B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RNF113B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	bA10G5.1; MGC26599; R113B; RING finger protein 113B; RNF113B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RNF113B

Image



Western blot analysis of extracts from COLO205/K562/Jurkat/293 cells, using RNF113B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.