




OSR2 Antibody

Product Code	CSB-PA976975
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2R0
Immunogen	Synthesized peptide derived from Internal of Human OSR2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total OSR2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FLJ90037; odd-skipped related 2 (Drosophila); OSR2; Protein odd-skipped-related 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OSR2

Image



Western blot analysis of extracts from Jurkat cells, using OSR2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.