




MNT Antibody

Product Code	CSB-PA107744
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99583
Immunogen	Synthesized peptide derived from internal of Human MNT.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MNT protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Max-binding protein MNT; Protein ROX; Myc antagonist MNT; MNT; ROX
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MNT

Image



Western blot analysis of extracts from A549 cells, using MNT antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.