




TEAD3 Antibody

Product Code	CSB-PA904251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99594
Immunogen	Synthesized peptide derived from internal of Human TEAD3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TEAD3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcriptional enhancer factor TEF-5; TEA domain family member 3; TEAD-3; DTEF-1; TEAD3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TEAD3

Image



Western blot analysis of extracts from Jurkat cells, using TEAD3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.