




TAOK3 Antibody

Product Code	CSB-PA935138
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2K8
Immunogen	Synthesized peptide derived from internal of Human TAOK3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TAOK3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase TAO3; Cutaneous T-cell lymphoma-associated antigen HD-CL-09; CTCL-associated antigen HD-CL-09; Dendritic cell-derived protein kinase; JNK/SAPK-inhibitory kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAOK3

Image



Western blot analysis of extracts from K562 cells, using TAOK3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.