




SIRT3 Antibody

Product Code	CSB-PA229200
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTG7
Immunogen	Synthesized peptide derived from C-terminal of Human SIRT3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SIRT3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NAD-dependent deacetylase sirtuin-3; mitochondrial; SIR2-like protein 3; hSIRT3; SIRT3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SIRT3

Image



Western blot analysis of extracts from COLO205 cells, using SIRT3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.