




ZNF446 Antibody

Product Code	CSB-PA235984
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWS9
Immunogen	Synthesized peptide derived from internal of Human ZNF446.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZNF446 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein 446; Zinc finger protein with KRAB and SCAN domains 20; ZNF446; ZKSCAN20;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF446

Image



Western blot analysis of extracts from HuvEc cells, using ZNF446 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.