




PLAGL1 Antibody

Product Code	CSB-PA256793
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM63
Immunogen	Synthesized peptide derived from internal of Human PLAGL1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PLAGL1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein PLAGL1; Pleiomorphic adenoma-like protein 1; Tumor suppressor ZAC; Lost on transformation 1; LOT-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLAGL1

Image



Western blot analysis of extracts from 293 cells, using PLAGL1 antibody.