




DGKB Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA290951 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y6T7 |
| Immunogen | Synthesized peptide derived from internal of Human DGKB. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous levels of total DGKB protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:3000 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | 90 kDa diacylglycerol kinase; DAG kinase beta; DAGK2; DGK-beta; DGKB |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | DGKB |

Image



Western blot analysis of extracts from HuvEC cells, using DGKB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.