




CD70 Antibody

Product Code	CSB-PA095181
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32970
Immunogen	Synthesized peptide derived from internal of Human CD70.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CD70 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CD70 antigen; CD27 ligand; CD27-L; Tumor necrosis factor ligand superfamily member 7; CD70
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD70

Image



Western blot analysis of extracts from HuvEC, K562, HepG2, Jurkat, HeLa and COLO205 cells, using CD70 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.