




FGF13 Antibody

Product Code	CSB-PA151977
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92913
Immunogen	Synthesized peptide derived from internal of Human FGF13.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total FGF13 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Fibroblast growth factor 13; FGF-13; Fibroblast growth factor homologous factor 2; FHF-2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FGF13

Image



Western blot analysis of extracts from Jurkat cells and HepG2 cells, using FGF13 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.