




HSPBAP1 Antibody

Product Code	CSB-PA920123
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EW2
Immunogen	Synthesized peptide derived from internal of Human HBAP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total HBAP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HSPB1-associated protein 1; 27 KdA heat shock protein-associated protein 1; Protein associated with small stress protein 1; HSPBAP1; PASS1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSPBAP1

Image



Western blot analysis of extracts from HeLa cells and 293 cells, using HBAP1 antibody.