




IL11RA Antibody

Product Code	CSB-PA043792
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14626
Immunogen	Synthesized peptide derived from C-terminal of Human IL11RA.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IL11RA protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Interleukin-11 receptor subunit alpha; IL-11 receptor subunit alpha; IL-11R subunit alpha; IL-11R-alpha; IL-11RA
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL11RA

Image



Western blot analysis of extracts from HeLa cells, using IL11RA antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.