




IL17RC Antibody

Product Code	CSB-PA185992
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NAC3
Immunogen	Synthesized peptide derived from internal of Human IL17RC.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IL17RC protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Interleukin-17 receptor C; IL-17 receptor C; IL-17RC; Interleukin-17 receptor-like protein; IL-17RL
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL17RC

Image



Western blot analysis of extracts from HeLa, HepG2 and HuvEc cells, using IL17RC antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.