




IFIT5 Antibody

Product Code	CSB-PA984140
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13325
Immunogen	Synthesized peptide derived from internal of Human IFIT5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IFIT5 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	IFIT-5; interferon-induced protein with tetratricopeptide repeats 5; retinoic acid- and interferon-inducible 58 kDa protein; RI58;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFIT5

Image



Western blot analysis of extracts from Jurkat cells, using IFIT5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.