




IL18R1 Antibody

Product Code	CSB-PA051971
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13478
Immunogen	Synthesized peptide derived from internal of Human IL18R.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IL18R protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CD218a; CDw218a; IL-1Rrp; IL1 receptor-related protein; IL18R
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL18R1

Image



Western blot analysis of extracts from Jurkat cells, using IL18R antibody.