




IFI6 Antibody

Product Code	CSB-PA253751
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09912
Immunogen	Synthesized peptide derived from C-terminal of Human IFI6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IFI6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	41806; FAM14C; G1P3; Ifi-6-16; IFI6
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFI6

Image



Western blot analysis of extracts from K562 cells, using IFI6 antibody.