




MYLK Antibody

Product Code	CSB-PA059969
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15746
Immunogen	Synthesized peptide derived from internal of Human MYLK.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total MYLK protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KMLS; KRP; MLCK; MYLK; myosin light chain kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYLK

Image



Western blot analysis of extracts from HepG2 cells, using MYLK antibody.