




CSRNP2 Antibody

Product Code	CSB-PA679642
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H175
Immunogen	Synthesized peptide derived from internal of Human TAIP-12.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TAIP-12 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TGF-beta-induced apoptosis protein 12; TAIP-12; Protein FAM130A1; Protein FAM130A1; Protein FAM130A1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CSRNP2

Image



Western blot analysis of extracts from HuvEC, K562, HepG2 and A549 cells, using TAIP-12 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.