




TRAF4 Antibody

Product Code	CSB-PA788984
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUZ4
Immunogen	Synthesized peptide derived from internal of Human TRAF4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TRAF4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TNF receptor-associated factor 4; Cysteine-rich domain associated with RING and Traf domains protein 1; Malignant 62; RING finger protein 83; TRAF4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRAF4

Image



Western blot analysis of extracts from 3T3 cells, using TRAF4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.