




ANAPC5 Antibody

Product Code	CSB-PA137607
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJX4
Immunogen	Synthesized peptide derived from internal of Human ANAPC5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ANAPC5 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ANAPC5; Anaphase-promoting complex subunit 5; ANC5; Similar to anaphase-promoting complex subunit 5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANAPC5

Image



Western blot analysis of extracts from HeLa cells, Jurkat cells and COLO205 cells, using ANAPC5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.