




CDC14A Antibody

Product Code	CSB-PA154120
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNH5
Immunogen	Synthesized peptide derived from Internal of Human CC14A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CC14A protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CC14A; CDC10 (cell division cycle 10; <i>S. cerevisiae</i> ; homolog); cdc14
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDC14A

Image



Western blot analysis of extracts from HeLa cells, using CC14A antibody.