




CHFR Antibody

Product Code	CSB-PA972040
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EP1
Immunogen	Synthesized peptide derived from Internal of Human CHFR.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CHFR protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	checkpoint with forkhead and ring finger domains; Checkpoint with forkhead and RING finger domains protein; CHFR; E3 ubiquitin-protein ligase CHFR; FLJ10796
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHFR

Image



Western blot analysis of extracts from Jurkat cells, using CHFR antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.